

測試時間僅需5秒!
驗血量極小!
操作容易!



Blood Glucose Monitoring System*

- 
Biosensor Technology
- 
Easy 1 Step Calibration
- 
5 Seconds Reading
- 
500 Memories Time & Date
- 
Average 7/14/30 Days
- 
1.0 µl
- 
Automatic Alarm Reminder Function
- 
Easy 1 Step Calibration

電子感應檢測試紙



4 個簡單步驟

-  **Match the code number**
Make sure the code on the display matches the vial or foil package of each test biosensor.
-  **Insert the test biosensor**
Insert only the bottom end of the test biosensor into the test meter.
-  **Apply sample**
Touch and hold a drop of blood to the top of the test biosensor.
-  **Take accurate result**
Read test result and pull off the test biosensor without blood contamination.