

# Intellisense Blood Pressure Monitor - Model EF-36W6052



2-user mega memory



- Clinically Proven accuracy
- **Advanced Position Sensor** (with flash and buzzer), hence for more accurate measurement readings
- Blood Pressure Level Indicator
- Detects body movement
- 100 sets readings / user with date and time
- Displays average of last three readings taken within last 10 minutes

## Advanced Position Sensor

When using wrist blood pressure monitor, measuring at heart level is important to achieve accurate blood pressure measurements.

The **Advanced Position Sensor** in EF-36W6052 ensures measurement begins only when the wrist is at the correct position.

While the position is incorrect, the Positioning indicator lights in orange. The buzzer beeps when the wrist is in proper position, and the Positioning indicator turns to blue. The wrist cuff automatically starts to inflate and measurement starts.



## Guide for proper measurement posture

Keep the position of the heart mark on the wrist cuff at the center of the chest.



香港浸信會醫院  
Hong Kong Baptist Hospital



浸信會愛羣社會服務處  
BAPTIST OI KWAN SOCIAL SERVICE